



Bipolar Cables

Silicone Bipolar Cable



Shafts with 3.0 m, 5.0m tip

LC - 190-100 10cm
LC - 190-105 16.5 cm

Silicone European Bipolar Cable



Shafts with 3.0 m, 5.0m tip

LC - 190-110 10cm
LC - 190-115 16.5 cm

Silicone European Bipolar Cable (Martin)



Shafts with 3.0 m, 5.0m tip

LC - 190-120 10cm
LC - 190-125 16.5 cm



Bipolar Cables

Silicone Bipolar Cable



Shafts with 3.0 m, 5.0m tip

LC - 190-130 10cm
LC - 190-135 16.5 cm

Silicone European Bipolar Cable (Erbe)



Shafts with 3.0 m, 5.0m tip

LC - 190-140 10cm
LC - 190-145 16.5 cm

Silicone Bipolar Cable (Erbe)



Shafts with 3.0 m, 5.0m tip

LC - 190-150 10cm
LC - 190-155 16.5 cm



Bipolar Cables


Silicone European Bipolar Cable (Valleylab)



Shafts with 3.0 m, 5.0m tip

LC - 190-160 10cm
LC - 190-165 16.5 cm


Silicone Bipolar Cable (Valleylab)



Shafts with 3.0 m, 5.0m tip

LC - 190-170 10cm
LC - 190-175 16.5 cm

Silicone Bipolar Artery Sealer Cables



Shafts with 3.0 m tip

LC - 190-180 10cm



Monopolar Cables

Silicone Monopolar Cable



Shafts with 3 m tip, 10cm long

LC - 190-230

(Instrument Connector: Screw Type Pin Plug Generator Connector Standar Bovie 8 mm Type)

Silicone Monopolar Endoscopic Cable



Shafts with 3 m tip, 10cm long

LC - 190-240

(Instrument Connector: 4 mm Female Fitting Generator Connector Standar Bovie 8 mm Plug)



Shafts with 10cm long

LC - 190-250

(Instrument Connector: 4 mm Female Fitting Generator Connector Standar Bovie 4 mm Banana Plug)



Monopolar Cables

Silicone Diathermy Cable



Shafts with 10cm long

LC - 190-260

(Instrument Connector: 4 mm Male(European Diathermy) Generator Connector Standar Bovie 4 mm Banana Plug)



Shafts with 10cm long

LC - 190-270

(Instrument Connector: 4.8 mm Female (English Diathermy) Generator Connector Standar Bovie 4 mm Banana Plug)



Shafts with 10cm long

LC - 190-280

(Instrument Connector: 4.8 mm Female(English Diathermy) Generator Connector Standar Standar 4 mm Banana Plug)